

ZMM-200 ABI METALLURGICAL MICROSCOPE



Inverted Metallurgical Microscope, Reflect, Polarizing Optional

- LWD Plan No Cover Glass Objective PL10x,20x,40x,100x
- Double Layer Mechanical Stage180x150mm, Moving Range 15x15mm,
- Built-in Round Plate On Working Stage Rotatable For Polarizing View
- Simple Polarizing Set Polarizer & Analyzer Slide Optional
- 6V20W Halogen Kohler Illumination, Brightness Adjustable, With Field Diaphragm, Aperture Diaphragm

Specifications	ZMM-200 ABI
Head	Binocular
Eyepiece	WF10x/18mm
Objective	LWD Plan (No Cover Glass) PL 10x/0.25 LWD Plan (No Cover Glass) PLL20x/0.40 LWD Plan (No Cover Glass) PLL40x/0.60 LWD Plan (No Cover Glass) PLL100x/1.25 (S,O)
Filter	Frosted Filter Blue Filter Green Filter Yellow Filter
Focusing	Coaxial Coarse & Fine Focus, Tension Adjustable, With Up Stop, Fine Focusing Scale 0.002mm
Nosepiece	Quadruple, Ball Bearing Inner Locating
Working Stage	Double Layer Mechanical, Size 180x150mm, Moving Range 15x15mm, With Built-in Round Plate On Working Stage Rotatable For Polarizing View
Light Source	6V20W Halogen Kohler Illumination, Brightness Adjustable, With Field Diaphragm, Aperture Diaphragm

Optional Accessories

Specifications	ZMM-200 ABI
Eyepiece	WF16x/11mm
	Dividing 10x, 0.1mm/Div
Objective	Plan Long Working Distance (No Cover Glass) PL5x/0.12
	Plan Long Working Distance (No Cover Glass) PLL60x/0.75(Spring)
	Plan Long Working Distance (No Cover Glass) PLL100x/0.85(Spring)
	Achromatic (No Cover Glass) 60x/0.85
Nosepiece	Quintuple Pall Pearing Inner Locating
Nosepiece	Quintuple, Ball Bearing Inner Locating
Simple Polarizing	Polarizer Slide, 90°Adjustable
	Analyzer Slide
CCD Adapter	0.4x
	0.5x
	1.0x
	0.5x With Dividing 0.1mm/Div
Photo Adapter	2.5x/4x Change Over Photograph Attachment With 10x Photo Eyepiece
	4x Focusing Photograph Attachment
	MD Adapter
	PK Adapter
	Digital Camera Adapter For Canon EF digital SLR Camera Adapter