

zygon

www.zygon.ca

info@zygon.ca



ZRM-100 AAA

ROLLER MIXER



Roller Mixer

Roller Mixer is primarily utilized for gentle mixing and maintaining biological samples in suspension. It is well-suited for various applications, such as preventing blood coagulation, latex diagnosis, immunoprecipitation, and more.

Specifications

ZRM-100 AAA

Speed range	20~99rpm
Operation mode	Rolling & Swinging
Timing range	1min~99h59min/∞
Forwards/Reverse motion	Yes
Timer for Forwards/Reverse	30sec~99min59sec
Amplitude	25mm
Roller quantity	6
Roller length	280mm
Max. load	4kg
Power	15W
Voltage	AC100~240V, 50/60Hz
External Size (W×D×H)(mm)	480x250x100
Shipping Size (W×D×H)(mm)	580X370x210
N.W./G.W. (kg)	3.8/4.7