

zygon

www.zygon.ca

| info@zygon.ca



ZRS-100 AAA ROCKING SHAKER



Rocking Shaker

Rocking shakers are suitable for mixing applications in wide fields including gel staining and destaining, washes, hybridisations, immune precipitation and Western blots etc.

- Double LCDs for independent display of time and speed.
- Different platforms and fixtures for selection.
- Smooth see-saw rocking motion with fixed angle of 9°.
- Wide range of platforms for use with a variety of vessels.

Specifications

ZRS-100 AAA

Rocking Motion	See-saw rocking
Angle	9°
Max. Load Capacity (With Attachment)	10 kg
Motor Type	DC Motor
Speed Range	10~70rpm
Speed And Timer Display	LCD
Time Range	1~1199 min
Operation Mode	Continuous/ Timed Operation
Working Environment	5-40°C, 80%RH
Motor Input	IP21
Motor Output	40W
Protection Type	24W
Consumption	50W

Power Supply	100~240V, 50/60Hz
Standard Accessory	1 pc petri dish platform
External Size (W × D × H)	360 × 410 × 200mm
Package Size (W × D × H)	540 × 480 × 340mm
Gross Weight	15 kg