

zygon

www.zygon.ca

| info@zygon.ca



ZTA-50 DB

TABLETOP AUTOCLAVE



Tabletop Autoclave

The unit disinfects and sterilizes the objects with saturated steam, by the means of pulsely exhausting the cold air to clear the efforts from which to the temperature, and make the objects completely drying by the way of vacuum-exhausting and jacket baking. By using the physical characteristics of the pressured saturated steam will releasing plenty of hot while condensing, the unit will keep the sterilizing objects wet and at high temperature in a given time, which spoils the virus and the spores together, and good sterilizing efforts will be received.

- The sterilizing procedure is automatically controlled, it is not necessary to be supervised, and used conveniently.
- At least 3 or more times pulse-vacuum exhausting, which is applied to wrapped or unwrapped objects, hollow equipment.
- Drying function.
- With inner printer, which records the full procedure of the sterilization.
- With Bowie&Dick testing function, which is for the steam penetrating testing.
- Self-testing malfunction, easy to maintain.

Specifications

ZTA-50 DB

Chamber Volume	50L
Max working temperature	134°C±1°C
Max working pressure	0.22mpa
Adjusted sterilizing temperature range	55°C~134°C
Adjusted timer range	1~99mins
Adjusted drying time range	1~99mins

Vacuum Limit	-0.08mpa
Hot Average	$\leq 2^{\circ}\text{C}$
Power Supply	AC220V/50HZ
Dimensions	770 × 600 × 540 mm
Transport Dimensions	880 × 730 × 700 mm
Gross/Net Weight	112/93 kg